

Expt. No. 3.

Integrated Optics  
Angle of Prism

Page No. 5

Apparatus Available:-

1. Spectrometer
2. Spirit level
3. Magnifying glass
4. Glass prism
5. Sodium vapour lamp

Objective:-

To determine the apex angle of given prism using a spectrometer

Calculation:-



~~$$\text{Mean } \theta = \frac{119^\circ 51' + 119^\circ 49'}{2} = 119^\circ 50'$$~~

~~$$\text{Angle of prism} = \frac{\text{Mean } \theta}{2} =$$~~

$$\text{Mean } \theta = \frac{120^\circ 44' + 121^\circ 13'}{2} = 120^\circ 58'$$

$$\text{Angle of prism} = \frac{120^\circ 58'}{2} = 60^\circ 29'$$

Teacher's Signature : \_\_\_\_\_



Reading of Reflected Ray	Vernier A			Vernier B		
	MSR	VSR	TOTAL	MSR	VSR	TOTAL
Reflection from side (a) AB	279°30'	7'	279°37'	99°	7'	99°7'
Reflection from side (b) AC	40° (+360°)	21'	400°21'	220°	20'	220°20'
Difference between a & b			120°44'			121°13'


*[Signature]*

*[Signature]* 20/12/18



Result :-

The apex angle of the given equilateral prism is  $60^{\circ} 29'$



Teacher's Signature : \_\_\_\_\_